



DuoSeal 1376 Kits

Major Repair Kit CAT. No. 1376K-06 for Model 1376

ID	Qty.	Description
Major Repair Kit 1376K-06		
1	1	Mechanical Shaft Seal Kit, 3/4" CAT. No. 1401D
2	1	Gasket, Intake Chamber Cover
3	1	Intake Screen
4	1	Thrust Washer
5	2	Woodruff Key
6	2	Woodruff Key
7	2	Gasket, Oil Case
8	1	O-Ring
9	1	Pressure Release Valve
10	1	Ball Bearing
11	1	Exhaust Valve, Lower
12	1	Exhaust Valve, Upper
13	2	Vane, Small
14	3	Spring Holder
15	6	Spring for Vanes
16	2	Vane, Large
17	4	Valve Spring
18	4	Valve for Block
19	2	Spacer for Valve Block

Minor Repair Kit CAT. No. 1376K-05

Contains same as Major Kit 1376K-06 except vane ID No. 13 and 16.