



DuoSeal 1374 Spare Parts

Individual Spare Parts for Model 1374

ID	Qty.	CAT. No.	Description
A	1	41-1492	Pulley with Set Screw
B	1	61-8770	Large End Plate
C	1	41-1484	Shaft and Collar
D	1	41-1478	Large Rotor
E	1	61-8691	Large Ring
F	1	1393L	Intake Nipple 1-5/8" I.D. Hose
G	1	41-1501	Dust Cap
H	1	41-1482	Oil Case Assembly
I	1	41-1500	Small End Plate
J	1	41-1497	Exhaust Ring
K	1	41-1499	Exhaust Rotor
L	1	61-8689	Center Plate
M	1	61-8693	Intake Chamber

1	1	★1401F	Mechanical Shaft Seal Kit, 1"
	1	★1374A	V-Belt - 1374
		★1417P-20	Standard Exhaust Filter
		★1417H-01	Replacement Exhaust Filter Element
		★1407K-11	DuoSeal Vacuum Pump Oil -1 Qt.
		★1407K-15	DuoSeal Vacuum Pump Oil -1 Ga.
		★331080-5	5' x 1-5/8" I.D. Vacuum Hose
5	1	41-1477	Large Vane Metal (Requires 2)
20	1	41-1487	Small Vane Metal (Requires 2)
	1	1371F-01	Belt Guard & Base Kit
	1	41-2814	Front Cover/Belt Guard
22	1	41-1736	Gas Ballast Valve
	1	41-1010	Motor Pulley with Set Screw
	1	41-3753	Oil Case Window (Requires 2)
	1	726130	Oil Drain Valve
	1	41-1925	Motor 1-1/2 HP, 115/230V, 60Hz, 1Ph