



ID	Qty.	CAT. No.	Description
A	1	41-2191	Pump Pulley w/ Set Screw
B	1	41-2369	Intake Chamber Cover
C	1	61-8476C	End Plate
D	1	41-2374	Shaft
E	1	41-2378	Oil Case Assembly
F	1	1393J	Hose Adapter for 7/16" I.D. Hose
G	1	41-1045	Front Plate
H	1	41-1040	Rotor
I	1	41-1042	Inner Ring
1	1	★1401E	Mechanical Shaft Seal Kit, 5/8"
	1	★1401E-01	Lip Shaft Seal Kit, 5/8"
	1	★1399A	V-Belt

ID	Qty.	CAT. No.	Description
		★1417	Standard Exhaust Filter
		★1417L	Replacement Exhaust Filter Element (2 pk.)
		★1407K-11	DuoSeal Vacuum Pump Oil -1 Qt.
		★1407K-15	DuoSeal Vacuum Pump Oil -1 Ga.
		★331040-5	5' x 7/16" I.D. Vacuum Hose

ID	Qty.	CAT. No.	Description
10	1	41-1059	Vane
	1	61-8901	Belt Guard
	1	61-8916	Pump Base
	1	41-2377	Motor Pulley with Set Screw
17	1	726130	Oil Drain Valve
	1	41-3753	Oil Case Window
	1	61-8622	1/3 HP Motor, 115V, 60Hz, 1Ph

DuoSeal 1399 Spare Parts