


DuoSeal 1399 Repair Kits



Major Repair Kit CAT. No. 1399K-04 for Model 1399

ID	Qty.	Description	ID	Qty.	Description
Major Repair Kit 1399K-04			10	2	Vane
1	1	Mechanical Shaft Seal Kit, 5/8" CAT. No. 1401E	11	2	Vane Spring
2	1	Gasket, Intake Chamber Cover	12	1	Spring Holder for Vane
3	1	Woodruff Key	13	1	Exhaust Valve
4	1	Woodruff Key	14	1	Lock Washer
5	1	Gasket Oil Case	15	1	Stud for Valve Cover
6	1	Screw	16	1	Flat Washer Brass
7	1	Thrust Disk			
8	1	Retaining Ring	Minor Repair Kit CAT. No. 1399K-03		
9	1	Washer for Shaft Collar	Contains same as Major Kit 1399K-04 except vane ID No. 10.		